

SAFETY INSTRUCTIONS

DANGER: TO REDUCE RISKS OF ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

- For adult use only.
- Use only on a dry clean surface of the body over clothing by lightly pressing and moving across the body for approximately 60 seconds per region.
- Use the ELEEEELS only on the body's soft tissue as desired without producing pain or discomfort. Do not use on head or any hard or bony area of the body.
- Use only the applicators that provide the best desired effect.
- Bruising can occur regardless of control setting or pressure applied. Check treated areas frequently and immediately stop at first sign of pain or discomfort.
- Keep fingers, hair or any other body part away from the shaft and backside of the applicator as pinching may occur.
- Do not place any objects into the ventilation holes of the ELEEEELS.
- Do not immerse in water or allow water to enter the ventilation holes of the ELEEEELS.
- Do not drop or misuse the ELEEEELS.
- Only recharge with supplied 24VDC HYPERICE charger.
- Carefully examine device and battery before each use.
- Do not tamper or alter ELEEEELS in any way.
- Never leave the ELEEEELS operating or charging unattended.



DISCLAIMER

PLEASE DO NOT USE THE ELEEEELS, OR ANY VIBRATION DEVICE WITHOUT FIRST OBTAINING APPROVAL FROM YOUR DOCTOR, SPECIALIST OR PHYSIOTHERAPIST especially if you are in pregnancy, have diabetes with complications such as neuropathy or retinal damage,

wear of pace-makers, recent surgery, epilepsy or migraines, herniated disks, spondylolisthesis, spondylol-ysis, or spondylosis, recent joint replacements or IUD's, metal pins or plates or any concerns about your physical health. Frail individuals and children should be accompanied by an adult when using any vibration device.

WARRANTY

IF WITHIN ONE YEAR FROM THE DATE OF PURCHASE, THIS PRODUCT FAILS DUE TO A DEFECT IN MATERIAL OR WORKMANSHIP, ELEEEELS WILL REPAIR OR REPLACE THE PRODUCT, OR NECESSARY COMPONENTS, FREE OF CHARGE.

THIS WARRANTY EXCLUDES:

- (a) damage caused by accident, abuse, mishandling, or transport;
- (b) units subjected to unauthorized repair;
- (c) units not used in accordance with ELEEEELS care instructions;
- (d) damage exceeding the cost of the product;
- (e) deterioration of the delivered product resulting from abnormal storage and/or safeguarding conditions on the client's premises, and
- (f) failure to show the dated proof of purchase.

X1 Percussive
Massage
Device



Unit 5, 15/F, Global Gateway Tower,
63 Wing Hong Street, Cheung Sha Wan,
Kowloon, Hong Kong
www.eleeeels.com

Instructions Manual



We develop ELEEEELS Percussive massage device to relieve pain and promote blood circulation to accelerate recovery and prevent muscle damage

CORDLESS/RECHARGEABLE
24 V Lithium-ion Battery

LIGHTWEIGHT DESIGN
2.5lbs / 1.1kg

SPECIFICATIONS



- 1 Powerful Brushless High-Torque Motor
- 2 Removable Battery Release
- 3 Rechargeable Lithium-Ion Battery (3+ Hrs of Use Per Charge)
- 4 6 Speed Settings for Optimal Pressure
- 5 Cordless and Lightweight (2.5 Lbs)
"Includes AC wall charger"
- 6 ON / OFF Power Switch
- 7 Interchangeable Applicators (x4)
- 8 ELEEEELS universal adaptor included

Vibration Frequency:
Comfort Mode: 900rpm, 1200rpm, 1500rpm
Sports Mode: 1800rpm, 2400rpm, 3200rpm

Weight: 2.5 lbs / 1.1 kg
Charge: 110-240V 50/60Hz 28V, DC-1A

USAGE INSTRUCTIONS

PLEASE READ PRIOR TO USE

CHARGING

- Fully charge battery for up to six hours before first use either attached or detached from the device
- To charge the battery:
 - 1) Battery can be charged when it is attached or detached from the device
 - 2) Power off the battery switch, and then connect the DC end of the supplied 28V adapter to the battery's charging port to charge
 - 3) The green light on the battery will blink during the charging process
 - 4) The green light stop blinking and remain fully illuminated to show it is fully charged
 - 5) Battery may be recharged at any time and at any battery level

OPERATING

- Align and attach battery into base of unit - ensure it is securely locked in place
- Power device on by moving the power switch on the bottom of the battery to the ON position. Power the device off by moving the power switch to the OFF position
- With power switch in ON position
- The massage gun will stop operation after a 10-minute session

Sports Mode:

Press the power button once to turn device into level 1 of the sports mode, further press will turn device to level 2 and 3 respectively. The blue LED lights will indicate the corresponding power level
User can press and hold the power button for 1 second to pause the device in working status

Comfort Mode:

Press and hold the power button for 2 seconds, and release the button will allow the device enter the level 1 of the comfort mode. Further press the button will allow the device to enter level 2 and 3 of the comfort mode. Again, User can press and hold the power button for 1

second to pause the device in working status

- Massage desired body part at the desired speed by applying pressure as needed
- To remove the battery, press the battery release button with thumb. Firmly pull downwards on the battery with the other hand
- The circular light indicator on the battery:
 - Red color : 0-25%
 - Blue color : 25-75%
 - Green color : 75- 100%

MAINTENANCE, CLEANING & STORAGE

- Use a slightly damp towel to clean device and a dry with a soft cloth
- Do not clean the device by chemicals or alcohols
- For storage or travel, make sure the power switch is powered off
- Only use the ELEEEELS supplied AC Adaptor

WARNINGS & CAUTION

- Do not immerse in water. Keep away from liquids or heat sources. Keep ventilation ports free from dust and debris.
- Do not remove screws or attempt to disassemble.
- Do not operate continuously for more than one hour. Allow device to rest for 10 minutes before re-using.
- The device will be halted for 30 minutes if it is continually operated over 1 hour with Sports Mode.
- Unplug the charger from device battery after charging or prior to use.
- For adult use only. Do not use if injured. Consult your doctor before using this product.